Arie Nasal Vented Mask
Mascarilla nasal Arie con ventilación
Masque nasal à fuite Arie
Arie Belüftete Nasenmaske
Arie neusmasker met ventilatieopeningen
Maschera nasale ventilata Arie
Máscara nasal com ventilação Arie
Arie ネーザルベントマスク
قناع أرفي بفتحات تنفيس

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Arie Nasal Vented Mask User Instructions

Arie Nasal Vented Mask is intended to be used with positive airway pressure devices, such as CPAP or bi-level, operating at or above 4 cm H₂O. The mask is to be used on adult patients (>30kg) for whom positive airway pressure therapy has been prescribed. The mask is intended for single-patient multi-use in the home, hospital, or institutional environment.

NOTES

- This mask does not contain PVC, DEHP or phthalates.
- This product is not manufactured with natural rubber latex.
- If the patient has ANY reaction to any part of the mask, discontinue use.
- An exhalation port is built into this mask, so a separate exhalation port is not required.
- The mask includes a 22mm connector compliant with ISO 5356-1 for compatibility with CPAP/Bi-level circuits.

CAUTIONS

- U.S. Federal Law restricts this device for sale by or on the order of a physician.
- Do not expose the Arie Nasal Vented Mask to temperatures above 140°F (60°C).
- This mask is intended for single patient use only. It may be cleaned and used repeatedly on the same person but must not be used on multiple persons. Do not sterilize or disinfect.

WARNINGS

- This mask is not suitable for providing life support ventilation.
- Do not use the Arie Nasal Vented Mask with baseline pressure less than 4 cm H₂O.
- If oxygen is used with this CPAP mask, the oxygen flow must be turned off when the CPAP machine is not operating. Explanation of the warning: When the CPAP device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP machine enclosure. Oxygen accumulated in the CPAP machine enclosure will create a risk of fire.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate.
- No smoking or open flames such as candles may be used when oxygen is in use.
- Any unusual skin irritation, chest discomfort, shortness of breath, gastric distension, abdominal pain, belching or flatulence from ingested air, or severe headache experienced by the patient during or immediately after use should be immediately reported to a healthcare professional.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult a physician or dentist.
- Consult a physician if any of the following symptoms occur while using the mask: drying of the eyes, eye pain, eye infections, or blurred vision. Consult an ophthalmologist if symptoms persist.
- Do not overtighten the headgear straps. Signs of overtightening include: excessive redness, sores, or bulging skin around the edges of the mask. Loosen headgear straps to alleviate symptoms.
- Significant unintentional leak may occur without properly fitting the mask. Follow the MASK FITTING INSTRUCTIONS provided.
- At low CPAP pressures, the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- This mask should be used with the positive airway pressure device recommended by a physician or respiratory therapist. The mask should not be used unless the positive airway pressure device is turned on and operating properly. Do not block or try to seal the exhalation port (vent holes). The vent holes allow a continuous flow of air out of the mask. When the positive airway pressure device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask's vent holes.
- This mask is not intended for patients without a spontaneous respiratory drive. This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask. This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia.

GENERAL INFORMATION

Dispose of the mask and accessories in accordance with local and regional guidelines.

- This mask has a non-transferable three-month limited warranty for manufacturer's defects from the date of purchase by the consumer. If a mask used under normal conditions should fail, Sleepnet will replace the mask. For additional information about the Sleepnet Warranty, please visit <http://www.sleepnethealth.com>.

com/. Per EU MDR 2017/745, please report serious incidents to manufacturer and EU Member State competent authority.

- Mask and headgear have a service life of 1 year.
- Make sure the mask is clean and dry before storing it for any period of time. Keep it in a dry place, away from direct sunlight.

ASSEMBLY INSTRUCTIONS

1. Insert the top of the tube into the hole at the front of the mask frame (Figure 1).
2. Align the V-shape of the frame with the V-shape slot of the mask cushion. Insert the frame into the mask cushion opening from the bottom to the top. Snap the frame into the mask by pushing both together (Figure 2).
3. Run your fingers along the connection to ensure the frame is flush with the mask cushion (Figure 3).
4. Connect the headgear clips to the mask frame by hooking the clips to the frame bars on the sides of the mask. Ensure the clips snap in securely (Figure 4).

FITTING INSTRUCTIONS

1. Detach one headgear clip from the frame (Figure 5).
2. Place the mask in the correct position on your face with the bottom of the mask cushion sealed just above your upper lip. The top of the cushion should rest on the bridge of your nose. Then reattach the headgear clip (Figure 6).
3. Tighten the top and bottom straps on both sides of the headgear so that the mask creates a snug seal on your face. Continue adjusting the mask until you achieve a comfortable seal (Figure 7).

DISASSEMBLY INSTRUCTIONS

1. Unhook the headgear clips from each side of the mask frame (Figure 8).
2. To remove the frame, gently push the top from inside while holding the cushion with your other hand (Figure 9).
3. Remove the tube from the mask frame. Grip and pull from the cuff to avoid stretching or damage (Figure 10).

CLEANING AND MAINTENANCE

Mask:

- Mask cleaning is recommended on a daily basis.
- Disassemble the mask as per the disassembly instructions above.
- Soak the components in warm water and mild liquid detergent solution for at least 2 minutes, but no more than 5 minutes (Figure 11).
- Gently clean the components with a soft bristle brush for at least 2 minutes, paying particular attention to the vent holes (Figure 12).
- Thoroughly rinse the components with fresh running water for at least 2 minutes to remove any residues.
- Inspect the components to ensure that they are visually clean. If visible soil remains, repeat the above cleaning steps.
- Allow the components to air dry out of direct sunlight (Figure 13).
- Re-assemble the mask as per the assembly instructions above.

Headgear:

- Headgear cleaning is recommended on a weekly basis.
- Hand wash the headgear in mild soap and warm water. Allow the headgear to air dry.

WARNING: Frequency of cleaning, methods of cleaning or the use of cleaning agents, other than those specified above, or exceeding the service life of the components can have an adverse effect on the device and consequently the safety or the quality of the therapy.

NOTE: Washing your face prior to fitting the mask will remove excess facial oils and help extend the life of the mask.

BEFORE EACH USE

- The mask and headgear should be replaced at the end of their service life of one year.
- Inspect the mask daily or prior to each use. Replace the mask if any component has material damage, such as tears, cracks or peeling of the cushion.

TECHNICAL INFORMATION

- Resistance to air flow at 50 L/min: 0.43 cm H₂O. Resistance to air flow at 100 L/min: 1.98 cm H₂O.
- Therapy pressure range: 4 to 30 cm H₂O.
- DEAD SPACE VOLUME (mL approx.) Small: 60 mL, Medium: 75 mL, Large: 95 mL.
- A-weighted Sound Pressure Level at 1 m distance: 20.4 dBA; A-weighted Sound Power Level: 28.4 dBA.
- Operating Temperature: 41°F to 104°F (5°C to 40°C).
- Operating Humidity: 15% to 95% RH non-condensing.
- Storage and transport temperature: -4°F to 140°F (-20°C to +60°C).

Sleepnet Arie is a trademark of Sleepnet Corporation. Visit our website at www.sleepnethealth.com.

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Medical Product Services GmbH