



CPAP

# Arie<sup>®</sup> Full Face Mask

**Light as air, strong in care**

The Sleepnet Arie<sup>®</sup> Full Face Mask is a next-generation CPAP mask engineered to improve sleep therapy through precision design, patient comfort, and reliability. Sleepnet innovation combines low-profile geometry, ergonomic headgear, and skin-friendly features to create a mask that patients can wear comfortably night after night.

## Clear field of vision


Low-profile design makes it easier to see.

## AIR<sup>°</sup>sil Technology

Our innovative AIR<sup>°</sup>sil technology is a soft, silicone rubber material engineered to provide lightweight comfort and a secure fit. Crafted from medical-grade silicone, AIR<sup>°</sup>sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

## Lightweight Frame

Minimizes bulk and enhances freedom of movement.



## 360° Rotating Swivel


Designed to provide greater flexibility and enhance comfort, the 360-degree rotating swivel adjusts to the wearer's natural head movement, allowing them to move freely in any direction during sleep.

## Replaceable Cushion

Our replaceable cushions can be removed and cleaned, helping to maintain better hygiene and mask performance.

**Available for Resupply**

S M L



## What makes the Arie® Mask different?

**Comfort** - Crafted from medical-grade silicone, AIR°sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

**Low Profile** - The Arie® Full Face Mask is a low profile option designed to provide a clear field of vision.

**Easy-to-use** - Easily fit while maintaining a reliable seal and minimizing skin irritation.

**Resupply** - The Arie® Full Face Mask uses a replaceable cushion, which is readily available for resupply.

### AIR°sil Technology

Innovative design engineered for comfort.



At Sleepnet, comfort and care are never an afterthought, they're the foundation of every mask we design. From our trusted AIR°gel® technology to our new AIR°sil innovation, we're redefining what restful breathing feels like.

### JustFit Headgear

Soft, ergonomic, and made for all-night wear.



· Four points of adjustment



· Two points of connection

### Arie® Full Face Mask with Headgear Specifications:

Product Number	Description	Cushion Size	Mask Weight
50804	Arie® Full Face Vented Mask Small (With headgear, frame, and tube)	Small	71.8 g
50805	Arie® Full Face Vented Mask Medium (With headgear, frame, and tube)	Medium	79.6 g
50806	Arie® Full Face Vented Mask Large (With headgear, frame, and tube)	Large	90.7 g
50851	Arie® Full Face Vented Mask Kit (Headgear, frame, tube, and S/M/L cushions)	S,M,L	N/A
50212	Arie® Full Face Vented Frame and Tubing Assembly	One-size	N/A
50231	Arie® Full Face Cushion	Small	34.4 g
50232	Arie® Full Face Cushion	Medium	41.7 g
50233	Arie® Full Face Cushion	Large	53.2 g
55043	JustFit Full Face Headgear	One-size	N/A



Also by Sleepnet®  
**Arie® Nasal**

### Contact us today:

☎ 1-800-742-3646  
✉ info@sleepnetcorp.com  
🌐 www.sleepnetmasks.com



The following product is protected by US patent pending

**FDA**  
Cleared

Does not contain natural rubber latex

M90215 Version 0



Sleepnet