



CPAP

Arie® Nasal Mask

Light as air, strong in care

The Sleepnet Arie® Nasal Mask is a next-generation CPAP mask engineered to improve sleep therapy through precision design, patient comfort, and reliability. Sleepnet innovation combines low-profile geometry, ergonomic headgear, and skin-friendly features to create a mask that patients can wear comfortably night after night.

Clear field of vision

Low-profile design makes it easier to see.

AIR°sil Technology

Our innovative AIR°sil technology is a soft, silicone rubber material engineered to provide lightweight comfort and a secure fit. Crafted from medical-grade silicone, AIR°sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

Lightweight Frame

Minimizes bulk and enhances freedom of movement.

360° Rotating Swivel


Designed to provide greater flexibility and enhance comfort, the 360-degree rotating swivel adjusts to the wearer's natural head movement, allowing them to move freely in any direction during sleep.


Replaceable Cushion

Our replaceable cushions can be removed and cleaned, helping to maintain better hygiene and mask performance.

Available for Resupply

S M L





What makes the Arie® Mask different?

Comfort - Crafted from medical-grade silicone, AIR°sil cushion technology gently adapts to the contours of each face, creating a soft, secure seal that minimizes air leaks and enhances comfort.

Low Profile - The Arie® Nasal Mask is a low profile option designed to provide a clear field of vision.

Easy-to-use - Easily fit while maintaining a reliable seal and minimizing skin irritation.

Resupply - The Arie® Nasal Mask uses a replaceable cushion, which is readily available for resupply.

AIR°sil Technology

Innovative design engineered for comfort.



At Sleepnet, comfort and care are never an afterthought, they're the foundation of every mask we design. From our trusted AIR°gel® technology to our new AIR°sil innovation, we're redefining what restful breathing feels like.

JustFit Headgear

Soft, ergonomic, and made for all-night wear.



· Four points of adjustment



· Two points of connection




Arie® Nasal Mask with Headgear Specifications:

Product Number	Description	Cushion Size	Mask Weight
50149	Arie® Nasal Vented Mask Small (With headgear, frame, and tube)	Small	40.9 g
50150	Arie® Nasal Vented Mask Medium (With headgear, frame, and tube)	Medium	43.8 g
50151	Arie® Nasal Vented Mask Large (With headgear, frame, and tube)	Large	45.7 g
50194	Arie® Nasal Vented Mask Kit (Headgear, frame, tube, and S/M/L cushions)	S,M,L	N/A
50218	Arie® Nasal Vented Frame and Tubing Assembly	One-size	N/A
50225	Arie® Nasal Cushion	Small	10.4 g
50226	Arie® Nasal Cushion	Medium	12.8 g
50227	Arie® Nasal Cushion	Large	15.0 g
55061	JustFit Nasal Headgear	One-size	N/A



Also by Sleepnet®
Arie® Full Face

Contact us today:

-  1-800-742-3646
-  info@sleepnetcorp.com
-  www.sleepnetmasks.com



The following product is protected by US patent pending

FDA
Cleared

Does not contain natural rubber latex

M90216 Version 0



Sleepnet